# OPERATING INSTRUCTIONS

## 60 L Water Bath

### 45-6490

<b>ELE International</b> Chartmoor Road, Chartwell Business Park Leighton Buzzard, Bedfordshire, LU7 4WG England phone: +44 (0) 1525 249200 fax: +44 (0) 1525 249249 email: <a href="mailto:ele@eleint.co.uk">ele@eleint.co.uk</a> <a href="http://www.ele.com">http://www.ele.com</a> <small>ELE International, a division of Hach Lange Ltd.</small>	<b>Distributor:</b>	<b>ELE International</b> Soiltest Product Division PO Box 389, Loveland, CO 80539 USA phone: +1 (800) 323 1242 fax: +1 (970) 663 9781 email: <a href="mailto:soiltest@eleusa.com">soiltest@eleusa.com</a> <a href="http://www.eleusa.com">http://www.eleusa.com</a>
<i>In the interests of improving and updating its equipment, ELE reserves the right to alter specifications to equipment at any time</i> <b>ELE International 2006 ©</b>		

## Contents

Section		Page
1	Safety	3
1.1	Key to hazard information	3
2	Introduction	3
3	Specification	3
4	Installation	4
5	Electrical Power Supply	4
6	Operation	4
7	Maintenance	5
8	Trouble shooting	5
9	Spares	5
WEEE Directive		

## 1 Safety

This equipment has been tested and is safe to use providing that the proper safety precautions are observed:

***Do not use this equipment in any manner, other than as specified in this user manual; misuse may result in serious injury to personnel.***

***Do not attempt to operate the equipment with covers removed.***

***Only connect to the correct electrical supply. Equipment voltage setting is stated on the fuse cover of the appliance inlet module.***

***Refer to Installation section before installing machine.***

***Do not operate machine with wet hands.***

Please read this entire manual before unpacking, setting up or operating this equipment. Pay particular attention to all **DANGER** and **CAUTION** statements. Failure to do so could result in serious injury to the operator, or other personnel, or damage to the equipment.

### 1.1 Key to hazard information

**DANGER :**

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION :**

Indicates a potentially hazardous situation that may result in minor or moderate injury.

**Important Note:** Information that requires special emphasis.

**Note:** Information that supplements points in the main text.

## 2 Introduction

- 2.1 The 60 L water bath is suitable for thermally conditioning Marshall samples prior to the Marshall stability test for both BS598 and EN12697-34 standards.
- 2.2 The bath is equipped with an advanced and accurate water temperature control system with a mechanical water circulating unit which will maintain a constant temperature of 60 +/- .5°C throughout the bath.
- 2.3 The water bath is supplied with a lid to minimise heat loss and water evaporation during use and the bath can condition up to 29 Marshall samples at one time.
- 2.4 For a full understanding of the conditioning process of Marshall samples, consultation of the required standard is recommended.

## 3 Specification

Maximum sample capacity	29
Supply voltage	230 VAC, 50/60Hz
Power consumption	1500W
Fuse	8 Amp (Fast blow) X 2
Maximum water capacity	60 L
Temperature range (Ambient)	5 - 40°C
Humidity range (Ambient)	30% - 90%
Certification	CE

## 4 Installation

- 4.1 When unpacking make sure the water bath is not damaged or any parts are missing.
- 4.2 It is recommended that the bath is located on a bench top near the samples to be stability tested, to minimise any drop in sample temperature when removed from the bath.
- 4.3 Fill the bath with water to the required level above the sample height before making any electrical connections to mains power. **Note:** during filling make sure the water does not accidentally access the top ventilation slats on the water bath.

## 5 Electrical Power Supply

### *Electrical Safety*

**Danger :** Before removing any cover or performing maintenance repair and service, isolate from electrical supply by removing mains plug. Where mains supply is required during these activities, only competent persons should perform the work.

**Danger :** Installation and servicing must be carried out by a qualified person.

**Danger :** A good low impedance Protective Earth Ground connection is required for this product.

**Danger :** It is recommended that this equipment is connected to a mains supply using an RCCB "Residual Current Circuit Breaker" for extra electrical safety protection.

**Danger :** Check that the power supply is compatible with the requirements stated on the voltage label on the unit and connect in accordance with IEE regulations or to local requirements.

**Danger :** Do not use the bath with wet hands. Dry hands carefully before touching any part of the bath or controls.

Power Cable is coded as follows:-

Brown wire	L	Live or Power
Blue wire	N	Neutral
Green/Yellow wire	E	Earth or Ground

## 6 Operation

**Caution :** Always use suitable protective garments to protect the operator from burns and scalds during use.

- 6.1 Never switch on the bath without it firstly being filled with water to the correct depth.
- 6.2 Switch on the bath using the green mains switch which will then illuminate. At the same time the LED display will show the current water temperature. The temperature of the water will slowly increase until it reaches the set temperature (see section 6.6 to change required set temperature). The stirrer unit should also be functioning.
- 6.3 Once the samples have been carefully placed into the bath it is essential to place the lid in position to minimise heat loss and evaporation of water.
- 6.4 During the heating process it is important to verify the water temperature of the bath using a glass thermometer. It might be necessary to slightly offset the set temperature to compensate for small errors. The lid has a small hole in one corner for this purpose.
- 6.5 Follow the required standard for temperature settings and sample soaking times.

- 6.6 The factory pre-set temperature is normally set to 37.8°C. To change the set temperature of the bath press the button marked "P" on the controller and, using the up and down arrows, adjust the display to the required set figure. By holding the arrows for an extended period this will speed up the increments. Once set the controller will automatically revert to the normal real-time temperature display.
- 6.7 During the sample conditioning process it is important to check that the height of the water has not dropped due to evaporation of the water. Replace as necessary.
- 6.8 Once the soaking process is finished and the samples removed, the bath can be switched off.

## **7 Maintenance**

**Danger** : Isolate the bath from the mains power supply before attempting any cleaning processes.

**Caution** : Make sure the water temperature has dropped to a safe level prior to any maintenance.

- 7.1 After repeated use it is recommended that sample debris is removed and the inside of the bath cleaned by initially removing the mains power supply, then removing the lower tray to access the bottom of the bath.

**Note:** Take care not to damage the heating element and thermostat probe.

- 7.2 The outer surfaces should only be cleaned with a soft dry cloth.
- 7.3 If the water bath is not used for any length of time, it is recommended that the mains power supply is disconnected and all the water drained from the bath prior to storage.

## **8 Trouble shooting**

If the mains switch does not illuminate when switched on, check the electrical power supply and plug for errors. Check the 2 fuses are intact.

If the switch light comes on but the bath does not heat up, the element could be damaged or the controller is not functioning correctly (contact ELE service department).

## **9 Spares**

Please contact ELE service department for any spares requirements.

## DIRECTIVE ON WASTE ELECTRICAL & ELECTRONIC EQUIPMENT (WEEE)

Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/96/EC), European electrical equipment users must now return old or end-of life equipment to the Producer for disposal at no charge to the user.

**Note:** For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment for proper disposal.

**Important document. Retain with product records.**

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