



IF IT'S WORTH BUILDING, IT'S WORTH TESTING

This wishlist was generated on 18/03/2020, and contains the following Products:

418-160

Paqualab System 50 for Analysis of Drinking Water to EC Standard Categories 1 & 2.



Paqualab System 50 for Analysis of Drinking Water to EC Standard Categories 1 & 2.

Code: **418-160**

Product Group: **Paqualab Standard Systems**

The two Paqualab systems are designed so that all the necessary equipment to determine the physical, chemical and microbiological parameters can be stored in one easily carried case. Results are obtained within 16 hours. The numerical results have precision and reproducibility. The method is simple and economical, enabling regular analysis to be carried out. This technique is accepted worldwide and Paqualabs are used in over 50 countries. Reusable components can be sterilised in the field. Consumables are pre-sterilised and have a long shelf life. Sample deterioration through transport and handling is eliminated.

Portable drinking water quality testing system. Microbiological - Coliform and E.Coli. Physical - temperature, conductivity and turbidity. Chemical - pH, Cl, NH₃, NO₂ and NO₃ plus over 40 chemical tests available. AC/DC and rechargeable battery powered incubator. On-site testing for fast, accurate results. Long-term monitoring and emergency testing of drinking water to EC and WHO guidelines. Comprising of the ELE Paqualab incubator 25 and 50, pH/temperature/mV meter, conductivity/temperature meter, turbidity meter, 0-1000 FTU, ELE Paqualab photometer, and consumables for 200 tests on faecal or total coliforms, chlorine, ammonia, nitrates and nitrites. The Paqualab is a portable testing system for the key drinking water quality parameters. The Paqualab allows the laboratory to be carried to water so that accurate results are obtained sooner, even in remote areas. Laboratory standard results can be obtained by operators with the minimum of training. The equipment provided in each Paqualab has been specifically designed for field use. Aluminium replaces glass for many components for improved strength and lighter weight. The incubators have 2 pre-set temperatures that are usually set to 37°C and 44°C. The internal battery provides sufficient power for a full incubation cycle. Over 1000 Paqualabs are in use throughout the world for monitoring drinking water in remote areas, emergency situations, research and education. The Paqualab is a completely modular system. Standard kits are available as shown below. All items are also available individually to enable the user to build a system to meet their specific requirements. In addition, Paqualab Incubator systems contain all the hardware required for microbiological tests. Paqualab standard systems are based on the Incubator systems and also include

tests for other parameters as shown below.

Further Information

Overall dimensions 475 x 230 x 275 mm.

Specification

Standard system

| Item supplied | 50 (418-160) | 25 (418-150) |
|---|--------------|--------------|
| (420-035) Universal & Incubator 50 | Yes | - |
| (420-030) Universal Incubator 25 | - | Yes |
| (430-020) pH/Temp/mV meter | Yes | Yes |
| (513-160) Conductivity/TDS/temp meter | Yes | Yes |
| (430-260) Turbidity meter | Yes | Yes |
| (430-550) Photometer | Yes | Yes |
| (433-115/10) Chlorine reagent system, free combined and total pack of 250 | Yes | Yes |
| (433-102/10) Ammonia reagent system, pack of 250 tests | Yes | Yes |
| (433-166/10) Nitrate reagent system, pack of 200 tests | Yes | Yes |
| (433-168/10) Nitrite reagent system, 4 packs of 50 tests | Yes | Yes |

Product Sheet

www.ele.com
+44 (0) 01525 249 200



| | | |
|---|--------------|--------------|
| (422-010) Coliform starter pack. Sufficient for 200 tests | Yes | Yes |
| Incubator system | | |
| Item supplied | 50 (418-160) | 25 (418-150) |
| Large rigid carrying case | Yes | - |
| Standard rigid carrying case | - | Yes |
| Universal incubator | Yes | Yes |
| Standard incubator | Yes | - |
| Incubator linking lead | Yes | - |
| Filtration unit | Yes | Yes |
| Mains lead | Yes | Yes |
| Battery lead | Yes | Yes |
| Cigar lighter plug | Yes | Yes |
| Aluminium petri dishes | Yes | Yes |
| Plastic autoclavable bottles (x4) | Yes | Yes |
| Dropping pipettes (x2) | Yes | Yes |
| Forceps | Yes | Yes |
| Magnifying glass | Yes | Yes |


| | | |
|------------------------------|-----|-----|
| Thermometer (spirit) | Yes | Yes |
| Trim tool | Yes | Yes |
| Grease | Yes | Yes |
| Instruction manual in binder | Yes | Yes |

Spares/Consumables




Aluminium Petri Dishes Pack of 25.

Code: [422-560](#)



Coliform Starter Pack. Sufficient for 200 Tests

Code: [422-010](#)



Filtration Unit AI Construction Complete with Sampling Cup.

Code: [420-450](#)



Membrane Lauryl Sulphate Medium 38.1 Gram Pack. Sufficient for 200 Tests.

Code: [422-110](#)



Pad Dispenser.

Code: [422-515](#)



Sterile Filter Membranes and Absorbent Pads. Pack of 200.

Code: [422-508](#)

Alternatives



Paqualab System 25 for Analysis of Drinking Water to EC Standard Categories 1 & 2.

Code: [418-150](#)