# OPERATING INSTRUCTIONS

## Heater/Magnetic Stirrer Unit

**46-4825**

<p><b>ELE International</b> Chartmoor Road, Chartwell Business Park Leighton Buzzard, Bedfordshire, LU7 4WG England phone: +44 (0) 1525 249200 fax: +44 (0) 1525 249249 email: <a href="mailto:ele@eleint.co.uk">ele@eleint.co.uk</a> <a href="http://www.ele.com">http://www.ele.com</a> ELE International, a division of Hach Lange Ltd.</p>	<p><b>Distributor:</b></p>	<p><b>ELE International</b> Soiltest Product Division PO Box 389, Loveland, CO 80539 USA phone: +1 (800) 323 1242 fax: +1 (970) 663 9781 email: <a href="mailto:soiltest@eleusa.com">soiltest@eleusa.com</a> <a href="http://www.eleusa.com">http://www.eleusa.com</a></p>
<p><i>In the interests of improving and updating its equipment, ELE reserves the right to alter specifications to equipment at any time.</i> <b>ELE International 2015 ©</b></p>		

---

## Contents

	<b>Section</b>	<b>Page</b>
	<b>1 Safety</b>	<b>3</b>
1.1	Key to hazard warning statements	<b>3</b>
	<b>2 Introduction</b>	<b>3</b>
	<b>3 Specification</b>	<b>3</b>
	<b>4 Installation</b>	<b>3</b>
	<b>5 Electrical Power Supply</b>	<b>3</b>
	<b>6 Operation</b>	<b>4</b>
	<b>7 Maintenance</b>	<b>5</b>
	<b>8 Trouble Shooting</b>	<b>5</b>
	<b>9 Spares</b>	<b>5</b>
	WEEE Directive	

## 1 Safety

This equipment has been tested and is safe to use providing that the proper safety precautions are observed:

**Do not use this equipment in any manner other than as specified in this user manual; misuse may result in serious injury to personnel.**

**Do not attempt to operate the equipment with covers removed.**

**Only connect to the correct electrical supply. Equipment voltage setting is stated on the label on the appliance.**

**Refer to Installation section before installing machine.**

**Do not operate the unit with wet hands.**

Please read this entire manual before unpacking, setting up or operating this equipment. Pay particular attention to all **DANGER** and **CAUTION** statements. Failure to do so could result in serious injury to the operator, or other personnel, or damage to the equipment.

### 1.1 Key to hazard warning statements

**DANGER :**

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION :**

Indicates a potentially hazardous situation that may result in minor or moderate injury.

**Important Note :** Information that requires special emphasis.

**Note :** Information that supplements points in the main text.

## 2 Introduction

The unit enclosed is used to heat samples for the ring and ball test. The unit has a variable speed magnetic stirrer, also a built in regulating transformer provides fine temperature control for the liquids used.

## 3 Specification

Supply voltage	230 VAC, 50/60 Hz
Power consumption	3.5 Amp
Plate temperature range (ambient)	50 - 350°C
Stirrer speeds	100-1800 rpm
Protection rating	IP 32
Weight	2.4 kg

## 4 Installation

4.1 When unpacking make sure the Heater/Magnetic Stirrer Unit is not damaged and that all parts are present.

4.2 It is recommended that the Heater/Magnetic Stirrer Unit is located on a stable bench top with a suitable 13 amp electricity power supply nearby.

## 5 Electrical Power Supply

### *Electrical Safety*

**Danger :** Before removing any covers or performing maintenance repair and service, isolate from the electrical supply by removing the mains plug from the supply. Where mains supply is required during these activities, only competent persons should perform the work.

**Danger** : Installation and servicing must be carried out by a qualified person.

**Danger** : A good low impedance Protective Earth Ground connection is recommended for this product.


**Danger** : Check that the power supply is compatible with the requirements stated on the voltage label on the unit and connect in accordance with IEE regulations or to local requirements.

Power Cable is coded as follows:-

Brown wire	L	Live or Power
Blue wire	N	Neutral
Green/Yellow wire	E	Earth or Ground

## 6 Operation

- 6.1 **Caution** : Always use suitable protective garments and gloves to protect the operator from burns and scalds during use of this product.
- 6.2 Control panel and power switch.

- 6.3 Place the required liquid in the beaker and place the beaker onto the heater platform.
- Note** : To get the correct temperature control during the test it might be necessary to use a heat transfer material between the beaker and heater platform.
- 6.4 Turn all the controls to their minimum positions and then switch on the mains switch on the front panel. Place the small stirrer bar into the liquid in the beaker and then, using the stirrer control knob, adjust to a suitable speed for stirring the liquid.
- 6.5 To start the heating process, adjust the temperature using the temperature control knob to add the required amount of heat to the liquid in the beaker.
- 6.6 **CAUTION** : Make sure that the liquid being heated does not get too hot

6.7 To stop the process switch off the mains switch on the unit.

6.8 **CAUTION** : Always isolate unit from mains electricity supply when not in use.

6.9 **Important Note** : If the Heater/Magnetic Stirrer Unit is being used only for heating applications then it is essential to have the stirrer motor running but without the stirrer bar in place to assist air circulation inside the case of the unit.

## 7 Maintenance

7.1 **Danger** : Isolate the Heater/Magnetic Stirrer Unit from the mains power supply before attempting any cleaning or maintenance processes.

7.2 **Caution** : Make sure the temperature has dropped to a safe level prior to any maintenance.

7.3 Clean the surface of the unit only using a damp cloth and mild detergent. Use a modern citrus based solvent for bitumen stains.

**Note** : For health and safety reasons the use of other solvent types is not recommended.

## 8 Trouble Shooting

If when the mains switch is switched on and no lights illuminate, check the electrical power supply and plug for errors. Also check all fuses are the correct capacity and are not blown. If the problem persists then contact the ELE service department for advice.

## 9 Spares

Please contact ELE service department for any spares requirements.

## DIRECTIVE ON WASTE ELECTRICAL & ELECTRONIC EQUIPMENT (WEEE)

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