



IF IT'S WORTH BUILDING, IT'S WORTH TESTING

This wishlist was generated on 10/05/2017, and contains the following Products:

---

42-0410

**Flakiness Gauge.**

---

42-0600

**Set of Flakiness Sieves Comprising 1 Ea 4.9mm 7.2mm 10.2mm 14.4mm 19.7mm 26.3mm and 33.9mm**

---

42-0820

**Length Gauge**

---

42-0821

**3.1 Shape Index Caliper.**

---

81-0590

**Vernier Caliper (LCD) Range 0 to 200mm X 0.01mm**

---

42-0300

**Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0302

**Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0304

**Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0306

**Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0308

**Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0310

**Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0314

**Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0316

**Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0318

**Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.**

---

---

42-0320

**Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0322

**Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0324

**Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.**

---

42-0326

**Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.**

---

46-5295/01

**Semi-Automatic Penetrometer**

---

46-5360/10

**Penetration Needle**

---

46-5800

**Transfer Dish**

---

46-5861

**Penetration Tin**

---

46-5860

**Penetration Tin**

---

45-6550/01

**30 L Heating Bath with Led Display Cover**

---

46-4605

**BS/EN Ring and Ball Apparatus**

---

46-4825/01

**Magnetic Stirrer**

---

82-5272

**Thermometer**

---

82-5274

**Thermometer**

---

45-3815/01

## Rotatest 3000

---

46-3310/01

## Cleveland Flash Cup Apparatus

---



## Flakiness Gauge.

Code: [42-0410](#)

Product Group: [Flakiness Sieves and Gauges](#)

Constructed of heavy gauge sheet steel to the dimensions specified in BS 812.

## Standards

BS 812

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Alternatives



## Length Gauge

Code: [42-0820](#)



## Set of Flakiness Sieves Comprising 1 Ea 4.9mm 7.2mm 10.2mm 14.4mm 19.7mm 26.3mm and 33.9mm

Code: [42-0600](#)

Product Group: [Flakiness Sieves and Gauges](#)

Comprising of 1 each, 4.9 mm, 7.2 mm, 10.2 mm,  
14.4 mm, 19.7 mm, 26.3 mm and 33.9 mm.

## Standards

BS 813

## Further Information


**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Alternatives



## Length Gauge

Code: [42-0820](#)



## Length Gauge

Code: [42-0820](#)

Product Group: [Elongation Index](#)

Manufactured to the dimensions specified in BS 182.

## Standards

BS 812-105

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Accessories



### 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



### 300mm Dia BS Sieve 14.0mm Perforated Plate.

Code: [79-2575](#)





## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)



## 300mm Dia BS Sieve 28.0mm Perforated Plate.

Code: [79-2615](#)



## 300mm Dia BS Sieve 37.5mm Perforated Plate.

Code: [79-2640](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)





## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.

Code: [79-2020](#)



## Air Jet Sieving Machine

Code: [80-0300/09](#)




## ELE Sieve Shaker

Code: [80-0200/06](#)



## ELE Sieve Shaker

Code: [80-0200/01](#)



## Sample Tray - 1200 x 1160 x 50 mm

Code: [81-4230](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)



## Shovel - With flat blade

Code: [81-0240](#)

## Alternatives




## Flakiness Gauge.

Code: [42-0410](#)



## Set of Flakiness Sieves Comprising 1 Ea 4.9mm 7.2mm 10.2mm 14.4mm 19.7mm 26.3mm and 33.9mm

Code: [42-0600](#)



## 3.1 Shape Index Caliper.

Code: [42-0821](#)

Product Group: [Determination of the Shape Index](#)

Designed to measure the ratio of length to width of individual aggregate particles.

## Standards

EN 933-4

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Accessories



### 200 mm Diameter Lid

Code: [79-0010](#)



### 200 mm Diameter Receiver

Code: [79-0020](#)



### 200mm Dia BS Sieve 1.00mm Stainless Steel Mesh.

Code: [79-0230](#)



## 200mm Dia BS Sieve 125Mic Stainless Steel Mesh.

Code: [79-0110](#)



## 200mm Dia BS Sieve 2.00mm Stainless Steel Mesh.

Code: [79-0270](#)



## 200mm Dia BS Sieve 250Mic Stainless Steel Mesh.

Code: [79-0150](#)



## 200mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-1500](#)



## 200mm Dia BS Sieve 500Mic Stainless Steel Mesh.

Code: [79-0190](#)

---



## 200mm Dia BS Sieve 63 Mic Stainless Steel Mesh.

Code: [79-0070](#)

---



## 200mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-1540](#)

---



## 300mm Dia BS Sieve 125mm Perforated Plate.

Code: [79-2710](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)

---





## 300mm Dia Receiver.

Code: 79-2020



## Air Jet Sieving Machine

Code: 80-0300/09



## ELE Sieve Shaker

Code: 80-0200/01



## ELE Sieve Shaker

Code: 80-0200/06



## Sample Tray - 406 x 406 x 50 mm

Code: 81-4030



## Sieve Brush Nylon

Code: [79-7210](#)



## Vernier Caliper (LCD) Range 0 to 200mm X 0.01mm

Code: [81-0590](#)

Product Group: [Determination of the Shape Index,](#)  
[Miscellaneous Laboratory Hardware](#)

Dual measuring range highly accurate vernier caliper 0 to 200mm x 0.01mm (8 inch x 0.001 inch) to comply with bs 887 and bs 6365. Complete with liquid crystal display.

## Standards

EN 933-4

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Accessories



### Vibrating Hammer


Code: [47-0455/01](#)

## Alternatives



### Vibrating Hammer

Code: [47-0455/01](#)



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)

Product Group: [Grid Sieves](#)

2.5mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 5 mm Retained Size - 4 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---



## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)





## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)

Product Group: [Grid Sieves](#)

3.15mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 6.3 mm Retained Size - 5 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---



## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)





## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)

Product Group: [Grid Sieves](#)

4mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 8 mm Retained Size - 6.3 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)

---



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)

---



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)

---



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)

---



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)

Product Group: [Grid Sieves](#)

5mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 10 mm Retained Size - 8 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0308](#)

Product Group: [Grid Sieves](#)

6.3mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 12.5 mm Retained Size - 10mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0306](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)

Product Group: [Grid Sieves](#)

8mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 16 mm Retained Size - 12.5 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)

---



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)

---



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)

---



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)

---



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)

---



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)

---



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)

---



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0306](#)

---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0308](#)



## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)

Product Group: [Grid Sieves](#)

10mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 20 mm Retained Size - 16 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0318](#)



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0316](#)

Product Group: [Grid Sieves](#)

12.5mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 25 mm Retained Size - 20 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)

---



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)

---



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)

---



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: 79-2010



## 300mm Dia Receiver.


Code: 79-2020



## Sample Tray - 406 x 406 x 50 mm


Code: 81-4030

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: 42-0314



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.


Code: 42-0318



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: 42-0306


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: 42-0308

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: 42-0310



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)

Product Group: [Grid Sieves](#)

16mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 31.5 mm Retained Size - 25 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)

---



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)

---



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)

---



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)

---



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)

Product Group: [Grid Sieves](#)

20mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 40 mm Retained Size - 31.5 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)

Product Group: [Grid Sieves](#)

25mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 50 mm Retained Size - 40 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)

---



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)

---



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)

---



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)

Product Group: [Grid Sieves](#)

31.5mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 63 mm Retained Size - 50 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)

---



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)

---



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)

---



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---




## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)

Product Group: [Grid Sieves](#)

40mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 80 mm Retained Size - 63 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: [79-2555](#)



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: [79-2560](#)



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---





## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)



## 300mm Dia Receiver.


Code: [79-2020](#)



## Sample Tray - 406 x 406 x 50 mm


Code: [81-4030](#)

## Alternatives




## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0316](#)



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0302](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.


Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)



## Semi-Automatic Penetrometer

Code: [46-5295/01](#)

Product Group: [Penetration Test](#)

As 46-5290 but incorporates Digital Automatic Controller which releases the needle assembly. The time set is displayed by a bright, easy to read display.

## Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Accessories



### Constant Temperature Bath 21 - 56°C

Code: [46-5500/01](#)



### Penetration Needle

Code: [46-5360/10](#)



### Penetration Needle

Code: [46-5340/C](#)



## Penetration Tin

Code: [46-5861](#)



## Transfer Dish

Code: [46-5800](#)

## Spares/Consumables



## ASTM (Non-Certified) Penetration Needle

Code: [46-5290/10](#)





## Penetration Needle

Code: [46-5360/10](#)

Product Group: [Penetration Test](#)

Hardened steel.

## Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## Transfer Dish

Code: [46-5800](#)

Product Group: [Bitumen Penetration Test](#), [Asphalt Penetrometer](#), [Penetration Test](#)

For use with ELE Penetrometer 46-5295.

## Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49



## Penetration Tin

Code: [46-5861](#)

Product Group: [Penetration Test](#)

For penetrations below 200, approximately 55 mm diameter x 35 mm deep.

## Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## Penetration Tin

Code: [46-5860](#)

Product Group: [Penetration Test](#)

For penetrations between 200 and 350,  
approximately 70 mm diameter x 45 mm deep.

## Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## 30 L Heating Bath with Led Display Cover

Code: [45-6550/01](#)

Product Group: [Accessories for manual and automatic compaction](#)

### Standards

BS 2000-49, EN 1426, EN 13179-2, ASTM D5, AASHTO T49, IP 49

### Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

### Accessories



#### Breaking Head (Marshall)

Code: [45-6850](#)

### Alternatives



#### Breaking Head (Marshall)

Code: [45-6850](#)





## BS/EN Ring and Ball Apparatus

Code: [46-4605](#)

Product Group: [Ring and Ball Apparatus](#)

With 2 shouldered pattern rings, 2 ball-centring guides and 2 balls. The apparatus has a support frame and is retained in a heat resistant container.

## Standards

EN 1427, BS 2000-58, UNI 4161, UNE 7111, DIN 52011, NFT66-008, AASHTO T53, ASTM D36, BS 2000

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

## Accessories




### Magnetic Stirrer

Code: [46-4825/01](#)



### Thermometer

Code: [82-5274](#)



## Thermometer

Code: [82-5272](#)





## Magnetic Stirrer

Code: [46-4825/01](#)

Product Group: [Ring and Ball Apparatus](#)

Utilises a rotating magnetic field to induce variable speed stirring action. A built-in regulating transformer provides fine temperature control of liquid up to a maximum temperature of +150°C. The unit is suitable for use with the Ring and Ball Apparatus.

## Standards

EN 1427, BS 2000-58

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## Thermometer

Code: [82-5272](#)

Product Group: [Ring and Ball Apparatus, Non Hazardous Glass Thermometers](#)

Mercury-free with a range of -10.0°C to + 210.0°C and graduated every 1°C. Immersion up to 76 mm

## Standards

EN 1427, BS 2000-58

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## Thermometer

Code: [82-5274](#)

Product Group: [Ring and Ball Apparatus, Non Hazardous Glass Thermometers](#)

Mercury-free with a range of -10.0°C to + 210.0°C and graduated every 1°C. Total immersion possible.

## Standards

EN 1427, BS 2000-58

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***



## Rotatest 3000

Code: [45-3815/01](#)

Product Group: [Asphalt Centrifuge Extractors](#)

Centrifuge extractor used to determine the quantitative amount of bitumen in bituminous paving mixtures. Variable speed control from 0 to 3600 rpm through the front panel mounted control knob and brake control included for rapid deceleration.

Precision machined aluminium removable cover with integral cup for adding solvent, with sealed cast aluminium housing.

For 220 V AC, 50 - 60 Hz, 1ph.

## Standards

EN 12697-1, ASTM D2172, AASHTO T164

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

3000 g capacity. Supplied with 100 Filter Discs. For 220 V AC, 50 - 60 Hz, 1ph.

## Accessories



### Bowl - 2.5 litre

Code: [81-5020](#)



### Desiccator - Non-vacuum

Code: [82-2100](#)



## Evaporating Dish - 100 x 40 mm

Code: [82-1970](#)



## Filter Discs Pack of 100

Code: [45-3807](#)



## Hotplate

Code: [78-2950/01](#)



## Measuring Cylinder - Glass 100 ml

Code: [82-0380](#)



## Measuring Cylinder - Glass 1000 ml

Code: [82-0500](#)



## Sample Tray - 306 x 306 x 38 mm

Code: [81-4020](#)



## Silica Gel

Code: [82-7091](#)



## Volumetric Flask - 100 ml

Code: [82-1000](#)



## Volumetric Flask - 1000 ml

Code: [82-1060](#)

## Alternatives



## Muffle Furnace

Code: [83-4170/01](#)

# Product Sheet

[www.ele.com](http://www.ele.com)  
+44 (0) 01525 249 200



## Rotatest 1500

Code: [45-3810/01](#)



## Cleveland Flash Cup Apparatus

Code: 46-3310/01

Product Group: [Cleveland Flash Cup Apparatus](#)

The ELE Semi-Automatic Cleveland Open Cup flashpoint Tester is equipped with an electrically heated cup (with a variable control to set temperature rise rate), a button operated and electrically driven sweep arm and a test flame for use on natural gas. Bench mounted, the unit is incorporated in an easy to clean stainless steel case.

## Standards

BS 2000-36, EN ISO 2592, ASTM D92, IP 36

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

For 220 - 240 V AC, 50 - 60 Hz, 1 ph.

## Specification

Temperature Range	Ambient to 400°C
Power	500 W maximum
Ramp rate	Variable, manually controlled
Voltage	220/240 V, 50/60 Hz
Size (HxWxD)	33 x 31 x 29 cm
Weight	6.5 kg