

















IF IT'S WORTH BUILDING, IT'S WORTH TESTING

Contents



This wishlist was generated on 21/07/2017, and contains the following Products:

36-3280/01

ADR Touch 2000 BS EN Compression Machine with Digital Readout and Self **Centring Platens.**

36-4150/01

ADR-Auto V2.0 2000 BS Compression Machine

37-5120

EN 12390-3/4 Distance Piece - 50 mm Depth

37-5140

EN 12390-3/4 Distance Piece - 60 mm Depth

37-5180

EN 12390-3/4 Distance Piece - 100 mm Depth



ADR Touch 2000 BS EN Compression Machine with Digital Readout and Self Centring Platens.

Code: 36-3280/01

Product Group: BS EN Compression Machines

The ADR Touch range of 2000 kN and 3000 kN capacity compression machines has been designed to meet the need for reliable and consistent testing. The load frame is a welded steel fabrication carrying the ball-seated upper platen. Positively located on the loading ram, which is protected from debris by a flexible cover, the lower platen is marked for the centring of cube and cylinder specimens. Self-centring lower platens for cube location are supplied as standard on EN machines and are available as an optional extra on the standard machine. The two machines for cube testing to EN standards are assembled and aligned using a special compression frame stability tester.

The dimensions of the frame allow the testing of concrete cylinders up to 320 mm long x 160 mm diameter, 150 and 100 mm square cubes, and on EN/BS machines, 200 mm square cubes. Kerbs and flagstones may also be tested on ADR machines as well as 150 mm and 100 mm square section beams to ASTM C78 using the optional 100 kN flexural frames which are connected to the power pack.

Standards

EN 12350-6, NF P18-411

Further Information

Not available for sale in the USA

ADR Touch BS EN 2000 and 3000 kN compression machines are supplied complete with self-centring lower platen and safety gates fitted with interlock switches ready for testing 300×150 mm diameter cylinders.

Whilst delivering all the features and reputation of the established ADR Series with its extensive design history, the new and improved user interface provides a high quality platform for testing that enhances the performance of ELE's compression machines.

New, sophisticated electronics further the benefits of simplified operation, whilst delivering the highest levels of accuracy in testing concrete and cement/mortar samples, satisfying the needs of Quality Control Managers, Lab Managers and Technicians.

ACCURACY & SAVINGS

The new ADR Touch Series, with 145mm (5.7") high resolution QVGA touch screen interface and intuitive menu-driven operation, reduces the time taken to set up the machine and perform tests, reducing the time to train staff by up to 25%.



Page 4 of 14

- Up to 6 sample types can be set as favorites, enabling one-touch set up for repeat testing.
- Full, QWERTY touch pad for input of test data.

Most testing errors produce lower strength results. Noncompliant loading rates can generate errors in measured strength. The user interface includes real-time display of load vs. time, further ensuring accurate and consistent test results and providing "goodness of test" data to improve traceability in your QC operations.

TRACEABILITY & DATA QUALITY

The ADR range now provides improved data quality and traceability in due diligence cases - it is now possible to demonstrate traceability all the way from the machine/user to the accreditation body, increasing your reputation and peace of mind - all test results now come complete with the machine serial number attached.

- Enhanced USB data communications between PC and machine eliminating the need for download software.
- Two gigabytes of storage memory.
- Full customization of sample sizes stress calculations are automatically recalculated.

USER SAFETY

With full safety gates as standard, total systems diagnostics, ram run-out switches and overload warnings ensure the

safety of your employees and the reliability of your machine.

Specification

Testing Standards BS EN 12390-4, BS EN

12390-3, BS EN 12390-4, BS EN 12390-5, EN 12504-1, EN 1354, EN 1521, EN 3161, EN 1338, EN 772-6, EN 13286-41, BS 3892-3, BS 187, BS 6717, BS EN ISO 7500-1, ASTM E4

Capacity 2000 kN
Cubes (Concrete) Up to 200 mm
Cylinders (Concrete) Up to 160 x 320 mm

Blocks Via optional Platen Handling

System

Flexural Testing Via Flexural Frame

TFV and ACV Yes

Voltage Supply 220-240V 50/60 Hz 1Ph.

Frame Class (Stability)

Frame Type

Max vert. Clearance

Max hor. Clearance

Welded

340 mm

355 mm

Platen Sizes Lower 220mm sq, Upper 300mm dia

Max Ram Travel 50 mm

Accessories

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242



100kN Flexural Fitting Kit (ADR) Used for Connecting Flexural Frames to ADR Compression Machines.

Code: 37-6138



2000 kN Capacity Calibration Load Cell with Hand-Held Readout

Code: 37-8315



ADR Touch Head for Compression machines

Code: 37-4950/09



BS/EN Standard Rectangular Platen

Code: 37-4860



EN 12390-3/4 Distance Piece - 100 mm Depth

Code: 37-5180



EN 12390-3/4 Distance Piece - 20 mm Depth

Code: 37-5110

www.ele.com +44 (0)1525 249 200

= +1 (800) 323 1242



EN 12390-3/4 Distance Piece - 50 mm Depth

Code: 37-5120



EN 12390-3/4 Distance Piece - 60 mm Depth

Code: 37-5140



Extended Front Safety Gate for Use with 2000kN BS/EN Load Frames Fitted with Rectangular Platens.

Code: 37-4835



Impact Printer Rs232 Serial Connection. Supplied Complete with 1 Paper Roll.

Code: 37-4859/01

Spares/Consumables



Pressure transducer, 0-700bar, 0.05 to 10v DC Output

Code: 6014A0062

www.ele.com +44 (0)1525 249 200



ADR-Auto V2.0 2000 BS Compression **Machine**

Code: 36-4150/01

Product Group: BS EN Compression Machines

The ADR-Auto V2.0 2000 BS EN is supplied complete with self-centring lower platen and safety gates fitted with interlock switches ready for testing 300 x 150 mm diameter cylinders. When used for cube testing Distance pieces (EN) of the appropriate size must be ordered separately.

Standards

EN 12350-6, NF P18-411

Further Information

Not available for sale in the USA

Whilst delivering all of the features and reputation of the established ADR-Auto V2.0 range with its 20 year design history, the new and improved user interface provides a high quality platform for testing that will enhance the performance of our compression machines. New, sophisticated electronics further the benefits of a closed-loop operation in testing concrete and cement/ mortar samples, satisfying the requirements of Quality Control Managers, Lab Managers and Technicians alike.

ACCURACY & SAVINGS

The NEW ADR-Auto V2.0 range, with 145mm (5.7") touch screen interface and intuitive menu driven operation reduces the time taken to set up the machine and perform tests, reducing the time to train staff by up to 25%. Most testing errors produce lower strength results. Noncompliant loading rates can generate errors in measured strength. With the latest closed-loop technology used in the ADR-Auto V2.0, the most effective cement ratios can be used with confidence.

TRACEABILITY & DATA QUALITY

The ADR-Auto V2.0 now provides improved data quality and traceability in due diligence cases - it is now possible to display traceability all the way from the machine/user to the accreditation body, increasing your reputation and peace of mind - all test results now come complete with the machine serial number attached. Prove your service quality - the ADR-Auto V2.0 now plots load vs time in real-time, proving full compliance to the testing standard. This data can be stored for later use or for detailed analysis of the concrete production.

USER SAFETY

Full safety gates as standard, total systems diagnostics, ram run-out switches and overload warnings ensure the safety of your employees and the reliability of your machine.

Specification

Testing Standards

BS EN 12390, BS EN 12390-3,

+44 (0)1525 249 200 **= +1 (800) 323 1242** Page 7 of 14 www.ele.com



BS EN 12390-4, BS EN 12390-5, EN 12504-1, EN 1354, EN 1521, EN 3161, EN 1338, EN 772-6, EN 13286-41, BS 3892-3, BS 187, BS 6717, BS EN ISO 7500-1, ASTM E4

Capacity 2000 kN
Cubes (Concrete) Up to 200 mm
Cylinders (Concrete) Up to 160 x 320 mm

Blocks Via optional Platen Handling

System

TFV and ACV Yes

Voltage Supply 220-240 V 50 Hz 1 Ph

Frame Type Welded
Max vert. Clearance 340 mm
Max hor. Clearance 355 mm

Platen Sizes Lower 220mm sq, Upper 300mm dia

Max Ram Travel 50 mm

Accessories



100 kN Flexural Fitting Kit ADR Auto for 220 - 240 V AC, 50 Hz, 1 ph

Code: 37-6135/01



100 kN Flexural Fitting Kit ADR Auto for 220 - 240 V AC, 60 Hz, 1 ph

Code: 37-6135/06



2000 kN Capacity Calibration Load Cell with Hand-Held Readout

Code: 37-8315

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 8 of 14



300mm Lower Brick Platen

Code: 37-4842



BS Block Platens & Platen Handling Assembly for 2000kN & 3000kN BS Frames.

Code: 37-4830



EN 12390-3/4 Distance Piece - 100 mm Depth

Code: 37-5180



EN 12390-3/4 Distance Piece - 20 mm Depth

Code: 37-5110



EN 12390-3/4 Distance Piece - 50 mm Depth

Code: 37-5120



EN 12390-3/4 Distance Piece - 60 mm Depth

Code: 37-5140

+44 (0)1525 249 200 www.ele.com

Page 9 of 14

Page 10 of 14



Flexible Hose Code: 5038B0027



Impact Printer Rs232 Serial Connection. Supplied Complete with 1 Paper Roll.

Code: 37-4859/01

Spares/Consumables



Pressure transducer, 0-700bar, 0.05 to 10v DC Output

Code: 6014A0062

Alternatives



ADR-Auto V2.0 2000 BS Compression Machine

Code: 36-4150/06



ADR-Auto V2.0 3000 BS EN Compression Machine

Code: 36-4165/01

www.ele.com +44 (0)1525 249 200 = +1 (800) 323 1242





ADR-Auto V2.0 3000 BS EN Compression Machine

Code: 36-4165/06

Page 12 of 14



EN 12390-3/4 Distance Piece - 50 mm Depth

Code: 37-5120

Product Group: Distance Pieces

Standards

EN 12390-3, BS 1881-115

Further Information

Not available for sale in the USA

Page 13 of 14



EN 12390-3/4 Distance Piece - 60 mm Depth

Code: 37-5140

Product Group: Distance Pieces

Standards

EN 12390-3, BS 1881-115

Further Information

Not available for sale in the USA



EN 12390-3/4 Distance Piece - 100 mm Depth

Code: 37-5180

Product Group: Distance Pieces

Standards

BS 1881-115, EN 12390-3

Further Information

Not available for sale in the USA