

















IF IT'S WORTH BUILDING, IT'S WORTH TESTING

Contents





This wishlist was generated on 05/03/2019, and contains the following Products:

24-8080/01

Automatic Soil Compactor BS CE Compliant

24-9200

BS Extension Collar

24-9202

BS Perforated Base Plate.

24-9204

BS Solid Base or Top Plate

24-9206

Cbr Cutting Collar with Area Ratio of 21.9% to Fit BS Moulds.

24-9208

C-Spanner for BS Cbr Mould. Two Required.

24-9210

Base Plate Tool for BS Mould

24-9212

Static Compaction Plug.

24-9220

Filter Papers 150mm Dia Equivalent to Whatman No 1. Box of 100.

24-9237

Spacer Disc 5-7/8" Diameter

25-1833/01

De-Aired Water Apparatus 15 Litre Capacity

25-3518/01

Digital Tritest 50 Load Frame

25-4047

50mm Triaxial Cell 1700Kpa with 5 Pressure/ Drainage Ports. Supplied with Two Valves.

25-4168

50mm Base Adapter for 50mm Cells.

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 2 of 36

Contents





25-4200

Piston Restraint Clamp for Ele Triaxial Cells

25-4290

Membrane Placing Tool 35mm/38mm.

25-4520

Valve No Volume Change 1/4 Inch BSp Fitted with 6 mm Connector and Integral Sealing Ring.

25-5050

Pressure Pad 38mm/1.5 Inch Diameter.

25-5061

Rubber Membrane 38mm/1.5 Inch Diameter. Pack of 10.

25-5081

Membrane Sealing Ring 1.5"

25-5100

Suction Membrane Device 38mm/1.5 Inch Diameter.

25-5120

Two-Way Split Former 38mm/1.5 Inch Diameter.

25-5430

Pressure Pad 50mm Diameter.

Rubber Membrane 50mm Diameter Pack of 10.

25-5461

Membrane Sealing Ring 50mm Diameter Pack of 10.

Membrane Placing Tool 50mm Diameter.

25-5480

Suction Membrane Device 50mm Diameter.

25-5500

Two-Way Split Former 50mm Diameter.

Contents



27-1798

DS7.2 California Bearing Ratio (CBR) Penetration Program for Windows 7, 32/64 bit

27-2798

DS7.3 CBR Penetation for Windows 10







Automatic Soil Compactor BS CE Compliant

Code: 24-8080/01

Product Group: BS CE Automatic Compactors

- With NEW CE Compliant guarding as an extra safety measure for users
- Five models meeting requirements of BS/EN and ASTM available
- Complies to BS 1377-4, BS 1924-2, EN 13286-2, EN 13286-47, ASTM D558, ASTM D698, ASTM D1557, AASHTO T99, AASHTO T134, AASHTO T180 standards
- Automatic re-setting of counter after completion of blow pattern

The time and effort required to prepare specimens for compaction studies and other test methods can often be costly and time-consuming. The use of an automatic, mechanical compactor will show considerable cost benefits over hand compaction methods.

These machines automatically compact specimens eliminating the laborious hand compaction method. The height and weight of the rammer are adjustable to suit test requirements.

An automatic blow pattern ensures optimum compaction for each layer of soil in compliance with the standards. The rammer travels across the mould and the table rotates the mould in equal steps on a base that is extremely stable. The number of blows per layer can be set at the beginning of the test.

Preset blow pattern ensures an even compaction, includes solid state controls for reliability and maintenance. An automatic digital counter resets to zero on the completion of a test. The selectable rammer has a drop height of 300-450mm, a selectable rammer weight of 2.5 or 4-5kg. A preset number of blows per layer can be set by thumb wheel control. A compaction rate of approximately 26 blows per minute. Accepts BS Standard Compaction and CBR moulds; also meets the requirements of BS 1377-4, BS 1924-2, EN 13286-2, EN 13286-47, ASTM D558, ASTM D698, ASTM D1557, AASHTO T99, AASHTO T134, AASHTO T180.

Specification Standards

BS 1377-4, BS 1924-2, EN 13286-2, EN 13286-47, ASTM D558, ASTM D698, ASTM D1557, AASHTO T99, AASHTO T134,

+44 (0)1525 249 200

www.ele.com +44 (0) 01525 249 200



AASHTO T180

Rammer BS/EN Circular faced, 50 mm dia,

selectable to 2.5 kg or 4.5 kg

weight

Rammer ASTM Circular faced, 2 in (50.8 mm)

dia, selectable to 5.5 lb (2.5 kg)

or 10 lb (4.5 kg) weight

Drop BS/EN Adjustable to 300 mm or 450 mm Drop ASTM


Adjustable to 12 in (305 mm)

or 18 in (457 mm)

Weight (kg) 165







BS Extension Collar

Code: 24-9200

Product Group: BS EN CBR mould accessories, Accessories for BS compaction and CBR


Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 7 of 36







BS Perforated Base Plate.

Code: 24-9202

Product Group: BS EN CBR mould accessories


To fit mould body

Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 8 of 36





BS Solid Base or Top Plate

Code: 24-9204


Product Group: BS EN CBR mould accessories, Accessories for BS compaction and CBR

Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 9 of 36





Cbr Cutting Collar with Area Ratio of 21.9% to Fit BS Moulds.

Code: 24-9206

Product Group: BS EN CBR mould accessories


To fit BS/EN mould body.

Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 10 of 36





C-Spanner for BS Cbr Mould. Two Required.

Code: 24-9208

Product Group: BS EN CBR mould accessories

To fit BS/EN CBR moulds and collars. Two required.


Standards

BS 1377-4, BS 1924-2, EN 13286-4

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 11 of 36







Base Plate Tool for BS Mould

Code: 24-9210

Product Group: BS EN CBR mould accessories

To fit into base plate to assist removal from mould body.


Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 12 of 36







Static Compaction Plug.

Code: 24-9212

Product Group: BS EN CBR mould accessories

Steel, 150 mm dia \times 51 mm thick. Complete with removable handle.

Standards

EN 13286-4, BS 1377-4, BS 1924-2

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 13 of 36





Filter Papers 150mm Dia Equivalent to Whatman No 1. Box of 100.

Code: 24-9220

Product Group: BS EN CBR mould accessories

Box of 100.

Standards

EN 13286-4, BS 1377-4, BS 1924-2







Spacer Disc 5-7/8" Diameter

Code: 24-9237

Product Group: Spacer Discs

5-7/8" diam. x 2" h. (147.6 x 50.8 mm).

Specification

Construction Finish Weights Steel; machined on all sides.

Plated.

Net 15.4 lbs. (7 kg).

www.ele.com +44 (0)1525 249 200





De-Aired Water Apparatus 15 Litre Capacity

Code: 25-1833/01

Product Group: De-Aired Water

Complies to BS 1377-1

This compact self-contained unit will de-air water quickly and efficiently down to levels of dissolved oxygen acceptable for geotechnical test methods. Air is removed from the water by a vacuum system, which continuously circulates the water in the tank. The unit is supplied with a clear water container, which will hold a maximum of 15 litres of water. Input and output lines are formed using standard 6 mm tube connectors.



Standards

BS 1377-1

Specification

Dimensions I x w x h (mm)

380 x 356 x 470

+44 (0)1525 249 200







Digital Tritest 50 Load Frame

Code: 25-3518/01

Product Group: Digital Tritest 50 Load Frame, CBR Load Frames, Triaxial Load Frames

- Microprocessor control
- Large on-board LED screen display
- Direct entry via a touch sensitive keyboard
- Rapid approach and return to datum of platen
- Fully variable speed, 0.00001 to 9.99999 mm/min
- Samples up to 100 mm diameter

This 50 kN capacity machine, designed primarily for triaxial testing of soil specimens up to 100 mm diameter x 200 mm long, comprises a rigid twin column construction with an integral fully variable microprocessor controlled drive unit and LCD display with a touch sensitive keyboard. The machine is normally bench mounted for ease of installation and operation.

The use of a microprocessor controlled drive system and keyboard entry provides the Digital Tritest 50 with a wide variety of features which include pause and speed reset during test, RS 232C, operator programming of speed and control functions, self test diagnostics and many other features.

A robustly constructed steel case houses the motor drive system with careful attention being given to the prevention of ingress of water or grit. All operating controls are mounted on the front panel of the machine, which is angled and recessed to prevent physical and environmental damage.

Standards

BS 1377-7, BS 1377-8, BS 1924-2, ASTM D1883, ASTM D2166, ASTM D2850, ASTM D4767, ASTM D7181, AASHTO T99, AASHTO T134, AASHTO T135, AASHTO T136, AASHTO T180, AASHTO T193, AASHTO T208, AASHTO T296, AASHTO T297

Further Information

Complete with RS 232C interface.

Specification

Capacity 50kN (11,200 lbf).

Speed Range English mode: 0.000001 to

0.399999 in/min.

Metric mode: 0.00001 to 9.99999 mm/min.

Rapid Approach Speed 1.0 in/min. (25 mm/min.).

www.ele.com # +44 (0)1525 249 200

www.ele.com +44 (0) 01525 249 200



Platen Travel 100mm (3.9"); limit switch

protection.

Platen Diameter 133 mm (5.2")

Vertical Clearance 910 mm, (36 .8") maximum; 305mm, (12") minimum.

Horizontal Clearance 364mm, (15.3").

Serial Interface RS232C; programmable baud

rate and protocol.

Overall Dimensions 500 x 500 x 1,470 mm, (19.7"

w. x 19.7" d. x 57.8" h).

Weight Net 100kg (220 lbs); Shpg.

136kg (300 lbs).

www.ele.com +44 (0)1525 249 200







50mm Triaxial Cell 1700Kpa with 5 Pressure/ Drainage Ports. Supplied with Two Valves.

Code: 25-4047

Product Group: Triaxial Cell, Triaxial Cells

- Working pressure up to 1700 kPa
- All round visibility
- Sample sizes 38 to 100 mm diameter
- Rapid assembly and dismantling
- Accepts a range of interchangeable submersible load transducers

This range of precision made triaxial cells has been designed to meet the requirements of the modern soils laboratory. The cells have been treated to minimise corrosion. Particular attention has been paid to the quality of finish between the piston and the head. Final assembly includes the fitting of an O-ring seal and the use of special lubricant to reduce friction to a minimum and eliminate water leakage.

The piston load capacity is designed to accept high horizontal forces which may be present during the final stages of a test. Each cell has five take-off positions drilled in the base for top drainage/back pressure, pore water pressure and bottom drainage. Two no-volume change valves and an anvil for strain gauge/transducer datum are supplied for fitting to the cell.

A feature of these cells is that they all accept a single diameter piston. The internal height is such that a range of submersible load transducers can be fitted without any modification. Each cell will accept a range of base adaptors and various accessories for testing a wide range of specimens

Standards

BS 1377, ASTM D2850, ASTM D4767, AASHTO T296, AASHTO T297

Specification

Cell size 50 mm

Max specimen size 50 x 100 mm

Working pressure 1700 kPa

Max piston load 45 kN

Vertical clearance required 380 mm

Horizontal clearance required 155 mm

Weight kg 4





50mm Base Adapter for 50mm Cells.

Code: 25-4168

Product Group: Specimen Base Adaptors, Base

Pedestals

50mm Base Adapter with Twin Pore Pressure Ports for 50mm Cells.

Standards

BS 1377, ASTM D2850, ASTM D4767, AASHTO T296, AASHTO T297







Piston Restraint Clamp for Ele Triaxial Cells

Code: 25-4200

Product Group: Aggressive Materials, General Cell

Accessories

For ELE Triaxial Cells Manufactured From July 1996.

Standards

BS 1377, ASTM D2850, ASTM D4767, AASHTO T296, AASHTO T297

Further Information

To lock the loading piston into position





Membrane Placing Tool 35mm/38mm.

Code: 25-4290

Product Group: Membrane Placing Tools, O-Ring

Placing Tools

Membrane Placing Tool For 35Mm/1.4 Inch And

38Mm/1.5 Inch Samples.

Standards

ASTM D2850, ASTM D4767, ASTM D5084, AASHTO T296, AASHTO T297, ASTM D7181

+44 (0)1525 249 200 **= +1 (800) 323 1242** Page 22 of 36 www.ele.com





Valve No Volume Change 1/4 Inch BSp Fitted with 6 mm Connector and Integral Sealing Ring.

Code: 25-4520

Product Group: Accessories for CU/CU effective stress, General Cell Accessories

Standards

ASTM D2850, ASTM D4767, ASTM D7181, AASHTO T296, AASHTO T297

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 23 of 36





Pressure Pad 38mm/1.5 Inch Diameter.

Code: 25-5050

Product Group: Pressure Pads

38 mm / 1.5 Inch diameter pressure pad. The pad shall have two empty ports and will come complete with two blanking plugs, plastic tubing and connector. Suitable For both drained and undrained tests.

Standards

ASTM D2850, ASTM D4767, ASTM D7181, AASHTO T296, AASHTO T297

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242

Page 24 of 36





Rubber Membrane 38mm/1.5 Inch Diameter. Pack of 10.

Code: 25-5061

Product Group: Rubber Membranes

Rubber Membrance For 38Mm/1.5 Inch Samples

(Pack Of 10).

Standards

AASHTO T296, AASHTO T297, ASTM D2850, ASTM D4767, ASTM D5084, ASTM D7181

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 25 of 36







Membrane Sealing Ring 1.5"

Code: 25-5081

Product Group: Sealing Rings

Membrane Sealing Ring 38Mm Diameter (Pack Of

10).

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 26 of 36







Suction Membrane Device 38mm/1.5 Inch Diameter.

Code: 25-5100

Product Group: Suction Membrane Devices, Suction

Membrane Device

Suction Membrane Device For 38Mm/1.5 Inch

Samples.





Two-Way Split Former 38mm/1.5 Inch Diameter.

Code: 25-5120

Product Group: Two-way Split Formers, Two-Way

Split Formers

Two-Way Split Former For 38Mm/1.5 Inch Samples.





Pressure Pad 50mm Diameter.

Code: 25-5430

Product Group: Specimen Caps, Pressure Pads

50 mm diameter pressure pad. The pad shall have two empty ports and will come complete with two blanking plugs, plastic tubing and connector. Suitable For both drained and undrained tests.

Standards

AASHTO T296, AASHTO T297, ASTM D2850, ASTM D4767, ASTM D5084, ASTM D7181





Rubber Membrane 50mm Diameter Pack of 10.

Code: 25-5441

Product Group: Rubber Membranes

50 mm diameter rubber membrane (Pack Of 10).

Standards

AASHTO T296, AASHTO T297, ASTM D2850, ASTM D4767, ASTM D5084, ASTM D7181

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 30 of 36





Page 31 of 36



Membrane Sealing Ring 50mm Diameter Pack of 10.

Code: 25-5461

Product Group: Membrane Sealing Rings, Sealing

Rings

50 mm diameter membrane sealing ring (Pack Of 10).

Standards

AASHTO T296, AASHTO T297, ASTM D2850, ASTM D4767, ASTM D5084

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242





Membrane Placing Tool 50mm Diameter.

Code: 25-5470

Product Group: Membrane Placing Tools, O-Ring

Placing Tools

50 mm Diameter Membrance Placing Tool.

Standards

AASHTO T296, AASHTO T297, ASTM D2850, ASTM D4767, ASTM D5084, ASTM D7181

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 32 of 36







Suction Membrane Device 50mm Diameter.

Code: 25-5480

Product Group: Suction Membrane Devices, Suction

Membrane Device

50 mm Diameter Suction Membrane Device.





Two-Way Split Former 50mm Diameter.

Code: 25-5500

Product Group: Two-way Split Formers, Two-Way

Split Formers

50 mm Diameter Two-Way Split Former.

www.ele.com +44 (0) 01525 249 200





DS7.2 California Bearing Ratio (CBR) Penetration Program for Windows 7, 32/64 bit

Code: 27-1798

Product Group: DS7.2 CBR Penetration Software, DS7.2 for Windows 7, Automatic Data Acquisition - CBR, DS7.2

This program provides the classical data for the CBR test. Load and penetration are monitored simultaneously. The program tabulates the results and provides a graphical presentation report.







DS7.3 CBR Penetation for Windows 10

Code: 27-2798

Product Group: DS7.3 for Windows 10, Automatic Data Acquisition - CBR, DS7.3

This program provides the classical data for the CBR test. Load and penetration are monitored simultaneously. The program tabulates the results and provides a graphical presentation report.