

















IF IT'S WORTH BUILDING, IT'S WORTH TESTING

Contents





This wishlist was generated on 12/04/2019, and contains the following Products:

24-4125/01

High Speed Stirrer Complete with Paddle Cup & Baffle. 220-240V 50/60Hz 1Ph.

24-4125/10

Mixing Paddle for 24-4125 Series High Speed Stirrer

24-4125/11

Soil Dispersion Cup & Baffle.

24-4145

Sodium Hexametaphosphate 500G.

24-4650

Soil Hydrometer Type 151H ASTM D422 Graduated 0.995 to 1.038 Specific Gravity.

24-4700

Hydrometer Sedimentation Cylinder

24-4865/01

Constant Temperature Bath 5 to 99.9 Deg.C X 0.1 Deg.C with Led Display and False Base Support.

79-5020

8in Dia ASTM Sieve 4.75mm Stainless Steel Mesh.

79-5070

8in Dia ASTM Sieve 2.0mm Stainless Steel Mesh.

79-5120

8in Dia ASTM Sieve 850 Mic Stainless Steel Mesh.

79-5160

8in Dia ASTM Sieve 425 Mic Stainless Steel Mesh.

79-5190

8in Dia ASTM Sieve 250 Mic Stainless Steel Mesh.

79-5240

8in Dia ASTM Sieve 106 Mic Stainless Steel Mesh.

79-5260

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 2 of 24

Contents



8in Dia ASTM Sieve 75 Mic Stainless Steel Mesh.

79-5400

8 Inch Dia Lid.

79-5510

8in Dia ASTM Sieve 9.50mm Stainless Steel Mesh.

79-5560

8in Dia ASTM Sieve 19.0mm Stainless Steel Mesh.

79-5580

8in Dia ASTM Sieve 25.0mm Stainless Steel Mesh.

79-5610

8in Dia ASTM Sieve 37.5mm Stainless Steel Mesh.

79-5630

8in Dia ASTM Sieve 50.0mm Stainless Steel Mesh.

79-5660

8in Dia ASTM Sieve 75.0mm Stainless Steel Mesh.







High Speed Stirrer Complete with Paddle Cup & Baffle. 220-240V 50/60Hz 1Ph.

Code: 24-4125/01

Product Group: Mechanical Analysis Stirrer

Compact bench top stirrer used for dispersing soil samples in water. Supplied complete with mixing paddle and dispersion cup with baffle. Beaker position maintained by column retainer.

Standards BS 1377, ASTM D422, AASHTO T88





Mixing Paddle for 24-4125 Series High Speed Stirrer

Code: 24-4125/10

Product Group: Mechanical Analysis Stirrer

Standards

BS 1377, ASTM D422, AASHTO T88





Soil Dispersion Cup & Baffle.

Code: 24-4125/11

Product Group: Mechanical Analysis Stirrer

Standards

BS 1377, ASTM D422, AASHTO T88







Sodium Hexametaphosphate 500G.

Code: 24-4145


Product Group: Sodium Hexametaphosphate,

Mechanical Analysis Stirrer

Standards

ASTM D422, AASHTO T88, ASTM E100





Soil Hydrometer Type 151H ASTM D422 Graduated 0.995 to 1.038 Specific Gravity.

Code: 24-4650

Product Group: ASTM Soil Hydrometers, Soil

Hydrometers

Graduated 0.995 to 1.038 g/ml.

Standards

ASTM D422, ASTM E100, AASHTO T88

Specification

Construction
Scale
Calibration

Bouyoucos; double scale. At 68°F.

Molded, seamless.

Weight Net 4 oz. (113 g).

+44 (0)1525 249 200







Hydrometer Sedimentation Cylinder

Code: 24-4700

Product Group: Hydrometer Sedimentation Cylinder,

Hydrometer Jar

Glass, 1000 ml capacity complete with rubber bung.

Standards BS 1337, ASTM D422, AASHTO T88

Further Information

1000Ml Glass







Constant Temperature Bath 5 to 99.9 Deg.C X 0.1 Deg.C with Led Display and **False Base Support.**

Code: 24-4865/01

Product Group: Constant Temperature Bath

This compact, bench-top stirrer is used for dispersing soil samples in water for hydrometer analysis. The stirrer is supplied complete with Mixing Paddle and Dispersion Cup. Weight 4 kg.

Further Information

Specially designed for the sedimentation testing of soils and other fine grained material, the bath is supplied with a false bottom to assist in circulation of the bath liquid. Will accommodate six Sedimentation Cylinders.

Specification

Dimensions (I x w x h) 535 x 210 x 610 mm external Capacity Holds up to 6 sedimentation

cylinders

Construction Stainless steel with toughened

glass front

Heater/thermostat/circulation Temperature Control

> with digital controller unit 5 to 99.9 Deg.C X 0.1 Deg.C

1500 watts Power

12 Weight kg

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 4.75mm Stainless Steel Mesh.

Code: 79-5020

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 4 Sieve Opening Standard 4.75 mm Sieve Opening Alternative (in.) 0.187

+44 (0)1525 249 200 www.ele.com **= +1 (800) 323 1242** Page 11 of 24







8in Dia ASTM Sieve 2.0mm Stainless Steel Mesh.

Code: 79-5070

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 10 2.00 mm Sieve Opening Standard Sieve Opening Alternative (in.) 0.0787

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 850 Mic Stainless Steel Mesh.

Code: 79-5120

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 20 Sieve Opening Standard 850 um Sieve Opening Alternative (in.) 0.331

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 425 Mic Stainless Steel Mesh.

Code: 79-5160

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 40 Sieve Opening Standard 425 um Sieve Opening Alternative (in.) 0.0165

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 250 Mic Stainless Steel Mesh.

Code: 79-5190

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Net 1 lb. (453 g). Weights

Sieve Size or Number No. 60 250 um Sieve Opening Standard Sieve Opening Alternative (in.) 0.0098

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 106 Mic Stainless Steel Mesh.

Code: 79-5240

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 140 106 um Sieve Opening Standard Sieve Opening Alternative (in.) 0.0041

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 75 Mic Stainless Steel Mesh.

Code: 79-5260

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 200 75 um Sieve Opening Standard Sieve Opening Alternative (in.) 0.0029

www.ele.com +44 (0)1525 249 200







8 Inch Dia Lid.

Code: 79-5400

Product Group: Sieve Lids and Receivers ASTM, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards ASTM E11







8in Dia ASTM Sieve 9.50mm Stainless Steel Mesh.

Code: 79-5510

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3/8" Sieve Opening Standard 9.5 mm Sieve Opening Alternative (in.) 0.375

+44 (0)1525 249 200 www.ele.com







8in Dia ASTM Sieve 19.0mm Stainless Steel Mesh.

Code: 79-5560

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3/4"
Sieve Opening Standard 19.0 mm
Sieve Opening Alternative (in.) 0.75

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 25.0mm Stainless Steel Mesh.

Code: 79-5580

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 1"

Sieve Opening Standard 25.0 mm

Sieve Opening Alternative (in.) 1

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 21 of 24







8in Dia ASTM Sieve 37.5mm Stainless Steel Mesh.

Code: 79-5610

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 1-1/2" Sieve Opening Standard 38.1 mm Sieve Opening Alternative (in.) 1.5

+44 (0)1525 249 200 www.ele.com **= +1 (800) 323 1242** Page 22 of 24







8in Dia ASTM Sieve 50.0mm Stainless Steel Mesh.

Code: 79-5630

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 2"

Sieve Opening Standard 50 mm

Sieve Opening Alternative (in.)

www.ele.com +44 (0)1525 249 200







8in Dia ASTM Sieve 75.0mm Stainless Steel Mesh.

Code: 79-5660

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

Standards

ASTM E11

Specification

Frames Full height; 2â^{III} (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3"

Sieve Opening Standard 75 mm

Sieve Opening Alternative (in.)

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242

Page 24 of 24