

















IF IT'S WORTH BUILDING, IT'S WORTH TESTING

# **Contents**





This wishlist was generated on 23/04/2019, and contains the following Products:

430-020

Digital pH/Temp/mV Meter with electrode, carrying case and pH 4 & 7 buffers

430-020/10

Spare pH electrode for use with 430-020 meter

46-3310/01

**Cleveland Flash Cup Apparatus** 

46-4605

**BS/EN Ring and Ball Apparatus** 

46-4825/01

**Electrical Hotplate with Integral Magnetic Stirrer** 

521-050

Dissolved Oxygen / Temp Meter. with Electrode and Carrying Case

78-2950/01

**Hotplate** 

79-5020

8in Dia ASTM Sieve 4.75mm Stainless Steel Mesh.

79-5060

8in Dia ASTM Sieve 2.36mm Stainless Steel Mesh.

79-5070

8in Dia ASTM Sieve 2.0mm Stainless Steel Mesh.

79-5100

8in Dia ASTM Sieve 1.18mm Stainless Steel Mesh.

79-5140

8in Dia ASTM Sieve 600 Mic Stainless Steel Mesh.

79-5160

8in Dia ASTM Sieve 425 Mic Stainless Steel Mesh.

79-5180

8in Dia ASTM Sieve 300 Mic Stainless Steel Mesh.

www.ele.com 👪 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 2 of 32

# **Contents**





79-5220

8in Dia ASTM Sieve 150 Mic Stainless Steel Mesh.

79-5260

8in Dia ASTM Sieve 75 Mic Stainless Steel Mesh.

79-5400

8 Inch Dia Lid.

79-5410

8in Dia Receiver.

79-5510

8in Dia ASTM Sieve 9.50mm Stainless Steel Mesh.

79-5530

8in Dia ASTM Sieve 12.5mm Stainless Steel Mesh.

79-5560

8in Dia ASTM Sieve 19.0mm Stainless Steel Mesh.

79-5580

8in Dia ASTM Sieve 25.0mm Stainless Steel Mesh.

79-5610

8in Dia ASTM Sieve 37.5mm Stainless Steel Mesh.

79-5630

8in Dia ASTM Sieve 50.0mm Stainless Steel Mesh.

79-5650

8in Dia ASTM Sieve 63.0mm Stainless Steel Mesh.

79-5660

8in Dia ASTM Sieve 75.0mm Stainless Steel Mesh.

79-7200

**Sieve Brush Brass and Nylon** 

82-5272

**Thermometer** 







# Digital pH/Temp/mV Meter with electrode, carrying case and pH 4 & 7 buffers

Code: 430-020

Product Group: Paqualab Test Meters

Robust waterproof case offering IP66/67 protection

Manual/automatic temperature compensation

Easy to use re-calibration function

The 430-020 pH meter is a three-in-one instrument that features a large easy to read, LCD display that indicates pH over the range of -2 to 16 pH with a resolution of 0.01 pH, mV over the range of -1000 to 1000 mV and temperature over the range of -39.9 to 149.9°C with a resolution of 0.1°C. The LCD display features both low battery indication and a user selectable backlight.

The pH readings are either manually or automatically temperature compensated over the range of 0 to 100°C. To automatically compensate, it is necessary to utilise a thermistor temperature probe. Each unit incorporates an auto-power off facility that automatically turns the instrument off after ten minutes, maximising battery life.

The 430-020 has an integrated rubber seal to ensure complete water resistance and helps to reduce the possiblity of damage in harsh environments. At the touch of a button, the instrument will automatically re-calibrate (two-point autocal) itself when used in conjunction with pH buffer solutions.

Each unit incorporates an easy to use BNC connector and Lumberg screw-locking type connector.

#### **Further Information**

#### Specification

Range Resolution Accuracy -2 to 16 pH 0.01 pH ± 0.02 pH ± 1000 mV 1 mV ±1 mV

-39.9 to 149.9 °C 0.1 °C ± 0.4 °C (-10 to 70 °C)

Battery 3 x 1.5 volt AAA

**Battery Life** Maximum 5 years (2500 hours) Combination electrode/ thermistor Sensor Type

Display 12 mm LCD **Dimensions** 32 x 71 x 141mm

+44 (0)1525 249 200

www.ele.com +44 (0) 01525 249 200



Weight

230 grams

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 5 of 32





Page 6 of 32



# Spare pH electrode for use with 430-020 meter

Code: 430-020/10

Product Group: Paqualab Test Meters

General purpose, epoxy bodied electrode with gel fill solution (max. temp. 60 deg. C.)

+44 (0)1525 249 200 **= +1 (800) 323 1242** www.ele.com







#### **Cleveland Flash Cup Apparatus**

Code: 46-3310/01

Product Group: Cleveland Flash Cup Apparatus

The ELE Semi-Automatic Cleveland Open Cup flashpoint Tester is equipped with an electrically heated cup (with a variable control to set temperature rise rate), a button operated and electrically driven sweep arm and a test flame for use on natural gas. Bench mounted, the unit is incorporated in an easy to clean stainless steel case.

#### **Standards**

BS 2000-36, EN ISO 2592, ASTM D92, IP 36

#### **Further Information**

\*\*NOT AVAILABLE FOR SALE IN THE USA\*\*

For 220 - 240 V AC, 50 - 60 Hz, 1 ph.

#### **Specification**

Temperature Range Ambient to 400°C
Power 500 W maximum
Ramp rate Variable, manually controlled
Voltage 220/240 V, 50/60 Hz
Size (HxWxD) 33 x 31 x 29 cm
Weight 6.5 kg

www.ele.com +44 (0)1525 249 200





# **BS/EN Ring and Ball Apparatus**

Code: 46-4605

Product Group: Ring and Ball Apparatus

With 2 shouldered pattern rings, 2 ball-centring guides and 2 balls. The apparatus has a support frame and is retained in a heat resistant container.

#### Standards

EN 1427, BS 2000-58, UNI 4161, UNE 7111, DIN 52011, NFT66-008, AASHTO T53, ASTM D36, BS 2000

#### **Further Information**

\*\*NOT AVAILABLE FOR SALE IN THE USA\*\*

www.ele.com +44 (0)1525 249 200







# **Electrical Hotplate with Integral Magnetic Stirrer**

Code: 46-4825/01

Product Group: Ring and Ball Apparatus

Circular-top analogue magnetic hotplate stirrer.

Speed controls for constant speed (even during load changes).

Anti-corrosive ceramic steel plate with excellent thermal conductivity.

Long life brushless motor.

High magnetic adhesion.

CE marked.

#### Standards

EN 1427, BS 2000-58

#### **Further Information**

\*\*NOT AVAILABLE FOR SALE IN THE USA\*\*

#### Specification

Temperature range Ambient-340 Deg.C

Stirring speed 0-1500 rpm

Max. Viscocity 40 cP/ mPas

Heating rate (1 litre water) 6K/min

Motor type Brushless

Stirring capacity 20L

Plate diameter 5.3"

Plate material Ceramic coated steel
Dimensions 11"L x 61/4"W x 3"H

Weight 6lbs

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 9 of 32





# Dissolved Oxygen / Temp Meter. with Electrode and Carrying Case

Code: 521-050

Product Group: Paqualab Test Meters

www.ele.com 🚆 +44 (0)1525 249 200 📕 +1 (800) 323 1242 Page 10 of 32







# Hotplate

Code: 78-2950/01 Product Group: Hotplate

#### **Further Information**

\*\*Not available for sale in the USA\*\*

# Specification

Construction Cast-iron top and base ventilated sheet-steel body

Nominal diameter 200 mm 500°C Maximum temperature Rated power 2000 W 5 kg

Weight







# 8in Dia ASTM Sieve 4.75mm Stainless Steel Mesh.

Code: 79-5020

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 4
Sieve Opening Standard 4.75 mm
Sieve Opening Alternative (in.) 0.187

www.ele.com # +44 (0)1525 249 200







# 8in Dia ASTM Sieve 2.36mm Stainless Steel Mesh.

Code: 79-5060

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 7
Sieve Opening Standard 2.36 mm
Sieve Opening Alternative (in.) 0.0937

www.ele.com +44 (0)1525 249 200







# 8in Dia ASTM Sieve 2.0mm Stainless Steel Mesh.

Code: 79-5070

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 10
Sieve Opening Standard 2.00 mm
Sieve Opening Alternative (in.) 0.0787

www.ele.com +44 (0)1525 249 200







# 8in Dia ASTM Sieve 1.18mm Stainless Steel Mesh.

Code: 79-5100

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

#### Standards

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 16 1.18 mm Sieve Opening Standard Sieve Opening Alternative (in.) 0.0469

www.ele.com +44 (0)1525 249 200 **=** +1 (800) 323 1242 Page 15 of 32







# 8in Dia ASTM Sieve 600 Mic Stainless Steel Mesh.

Code: 79-5140

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 30 Sieve Opening Standard 600 um Sieve Opening Alternative (in.) 0.0234

www.ele.com +44 (0)1525 249 200







# 8in Dia ASTM Sieve 425 Mic Stainless Steel Mesh.

Code: 79-5160

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 40 Sieve Opening Standard 425 um Sieve Opening Alternative (in.) 0.0165

www.ele.com +44 (0)1525 249 200







# **8in Dia ASTM Sieve 300 Mic Stainless** Steel Mesh.

Code: 79-5180

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

# Standards

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Net 1 lb. (453 g). Weights

Sieve Size or Number No. 50 300 um Sieve Opening Standard Sieve Opening Alternative (in.) 0.0117

www.ele.com +44 (0)1525 249 200







# 8in Dia ASTM Sieve 150 Mic Stainless Steel Mesh.

Code: 79-5220

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

# Standards

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 100 Sieve Opening Standard 150 um Sieve Opening Alternative (in.) 0.0059

www.ele.com +44 (0)1525 249 200







# 8in Dia ASTM Sieve 75 Mic Stainless Steel Mesh.

Code: 79-5260

Product Group: 8in Dia ASTM Brass Sieves, 8"
Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

**Diameter ASTM Sieves** 

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number No. 200 Sieve Opening Standard 75 um Sieve Opening Alternative (in.) 0.0029

www.ele.com +44 (0)1525 249 200







# 8 Inch Dia Lid.

Code: 79-5400

Product Group: Sieve Lids and Receivers ASTM, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards ASTM E11

+44 (0)1525 249 200







# 8in Dia Receiver.

Code: 79-5410

Product Group: Sieve Lids and Receivers ASTM, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

Standards ASTM E11







# 8in Dia ASTM Sieve 9.50mm Stainless Steel Mesh.

Code: 79-5510

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# Specification

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3/8"
Sieve Opening Standard 9.5 mm
Sieve Opening Alternative (in.) 0.375

www.ele.com +44 (0)1525 249 200







# **8in Dia ASTM Sieve 12.5mm Stainless** Steel Mesh.

Code: 79-5530

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

#### Standards

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 1/2" Sieve Opening Standard 12.5 mm Sieve Opening Alternative (in.) 0.5

+44 (0)1525 249 200 www.ele.com **=** +1 (800) 323 1242 Page 24 of 32







# **8in Dia ASTM Sieve 19.0mm Stainless** Steel Mesh.

Code: 79-5560

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

# Standards

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3/4" Sieve Opening Standard 19.0 mm Sieve Opening Alternative (in.) 0.75

+44 (0)1525 249 200 www.ele.com







# **8in Dia ASTM Sieve 25.0mm Stainless** Steel Mesh.

Code: 79-5580

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

#### Standards

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number

25.0 mm Sieve Opening Standard

Sieve Opening Alternative (in.)

+44 (0)1525 249 200 www.ele.com

**=** +1 (800) 323 1242

Page 26 of 32







# 8in Dia ASTM Sieve 37.5mm Stainless Steel Mesh.

Code: 79-5610

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 1-1/2"
Sieve Opening Standard 38.1 mm
Sieve Opening Alternative (in.) 1.5

www.ele.com # +44 (0)1525 249 200







# 8in Dia ASTM Sieve 50.0mm Stainless Steel Mesh.

Code: 79-5630

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 2"

Sieve Opening Standard 50 mm

Sieve Opening Alternative (in.) 2

www.ele.com +44 (0)1525 249 200 = +1 (







# **8in Dia ASTM Sieve 63.0mm Stainless** Steel Mesh.

Code: 79-5650

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8"

Diameter ASTM Sieves

#### Standards

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacquered brass

with S.S. mesh.

Show standard, alternate and Nameplates

Tyler Screen Scale equivalent

designations.

Net 1 lb. (453 g). Weights

Sieve Size or Number 2-1/2" Sieve Opening Standard 63 mm Sieve Opening Alternative (in.) 2.5

+44 (0)1525 249 200 www.ele.com **=** +1 (800) 323 1242 Page 29 of 32







# 8in Dia ASTM Sieve 75.0mm Stainless Steel Mesh.

Code: 79-5660

Product Group: 8in Dia ASTM Brass Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves, 8" Diameter ASTM Sieves

#### **Standards**

ASTM E11

# **Specification**

Frames Full height; 2â<sup>III</sup> (51 mm) depth

to mesh.

Construction Polished and lacguered brass

with S.S. mesh.

Nameplates Show standard, alternate and

Tyler Screen Scale equivalent

designations.

Weights Net 1 lb. (453 g).

Sieve Size or Number 3"

Sieve Opening Standard 75 mm

Sieve Opening Alternative (in.) 3

www.ele.com +44 (0)1525 249 200 = +1 (80







# **Sieve Brush Brass and Nylon**

Code: 79-7200

Product Group: Sieve Brushes

Double-ended, brass and nylon bristle

# Standards ASTM E11

#### **Further Information**

\*\*NOT AVAILABLE FOR SALE IN THE USA\*\*

www.ele.com +44 (0)1525 249 200

www.ele.com +44 (0) 01525 249 200





#### **Thermometer**

Code: 82-5272

Product Group: Non Hazardous Glass Thermometers,

Ring and Ball Apparatus

Mercury-free with a range of -10.0°C to + 210.0°C and graduated every 1°C. Immersion up to 76 mm

#### **Standards**

EN 1427, BS 2000-58

#### **Further Information**

\*\*NOT AVAILABLE FOR SALE IN THE USA\*\*