

## Grid Sieve 3.15mm Slot Width for EN Aggregate Flakiness Test.

Code: 42-0302

Product Group: Grid Sieves

3.15mm

## Standards

EN 933-3

## Further Information

**\*\*NOT AVAILABLE FOR SALE IN THE USA\*\***

Passing Size - 6.3 mm Retained Size - 5 mm

## Accessories



## 300mm Dia BS Sieve 10.0mm Perforated Plate.

Code: 79-2555



## 300mm Dia BS Sieve 11.2mm Perforated Plate.

Code: 79-2560



## 300mm Dia BS Sieve 12.5mm Perforated Plate

Code: [79-2565](#)

---



## 300mm Dia BS Sieve 16.0mm Perforated Plate.

Code: [79-2580](#)

---



## 300mm Dia BS Sieve 20.0mm Perforated Plate.

Code: [79-2595](#)

---



## 300mm Dia BS Sieve 25.0mm Perforated Plate.

Code: [79-2605](#)

---



## 300mm Dia BS Sieve 31.5mm Perforated Plate.

Code: [79-2630](#)

---



## 300mm Dia BS Sieve 4.00mm Perforated Plate.

Code: [79-2500](#)

---



## 300mm Dia BS Sieve 40.0mm Perforated Plate

Code: [79-2645](#)

---



## 300mm Dia BS Sieve 45.0mm Perforated Plate

Code: [79-2650](#)

---



## 300mm Dia BS Sieve 5.00mm Perforated Plate.

Code: [79-2515](#)

---



## 300mm Dia BS Sieve 5.60mm Perforated Plate.

Code: [79-2520](#)

---



## 300mm Dia BS Sieve 50.0mm Perforated Plate.

Code: [79-2655](#)

---



## 300mm Dia BS Sieve 6.30mm Perforated Plate.

Code: [79-2525](#)

---



## 300mm Dia BS Sieve 63.0mm Perforated Plate.

Code: [79-2670](#)

---



## 300mm Dia BS Sieve 8.00mm Perforated Plate.

Code: [79-2540](#)

---



## 300mm Dia BS Sieve 80mm Perforated Plate

Code: [79-2684](#)

---



## 300mm Dia BS Sieve 90.0mm Perforated Plate.

Code: [79-2690](#)

---



## 300mm Dia Lid.

Code: [79-2010](#)

---



## 300mm Dia Receiver.

Code: [79-2020](#)

---




## Sample Tray - 406 x 406 x 50 mm

Code: [81-4030](#)

---


## Alternatives



## Grid Sieve 10mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0314](#)


---



## Grid Sieve 12.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0316](#)


---



## Grid Sieve 16mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0318](#)


---



## Grid Sieve 2.5mm Slot Width for EN Aggregate Flakiness Test

Code: [42-0300](#)


---



## Grid Sieve 20mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0320](#)


---



## Grid Sieve 25mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0322](#)


---



## Grid Sieve 31.5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0324](#)


---



## Grid Sieve 40mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0326](#)


---



## Grid Sieve 4mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0304](#)


---



## Grid Sieve 5mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0306](#)


---



## Grid Sieve 6.3mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0308](#)

---



## Grid Sieve 8mm Slot Width for EN Aggregate Flakiness Test.

Code: [42-0310](#)