
### **Product Sheet**







#### **Shear Box Assembly 60mm Sq.**

Code: 26-2181

Product Group: Shear Box Assemblies, Shear Box Assemblies

- •Manufactured from corrosion resistant materials
- •Sample size 60 x 60mm x 25mm high meeting BS 1377
- •Incorporates PTFE faced box separation screws to reduce friction
- •Supplied complete with upper and lower porous plates and loading pad
- •Suitable for guick undrained or drained/residual tests

#### **Standards**

BS 1377, EN 1997-2, ASTM D3080, AASHTO T236

#### **Further Information**

All shearbox assemblies are designed to fit the carriage of 26-2114 series Digital Direct/Residual Shear Apparatus. They are supplied complete with three porous plates, one retaining plate and a loading pad. The three sizes supplied comply with the relevant requirements of BS 1377 and ASTM D3080. All assemblies can be used for guick shear tests or drained/residual shear tests.

Optional accessories are available including specimen cutters and extrusion tools.

#### **Specification**

Specimen area 60 x 60 mm

Specimen thickness 25 mm

Weight kg 2

Relevant standard BS 1377

#### **Accessories**



2.42" Diameter Shear Box Assembly (California Type)

Code: 26-2223

www.ele.com # +44 (0)1525 249 200 # +1 (800) 323 1242 Page 1 of 3

# **Product Sheet**







#### **Shear Box Assembly 100mm Sq.**

Code: 26-2197



### Specimen Cutter 60mm Sq.

Code: 26-2185



### Specimen Extrusion Tool 60 mm sq.

Code: 26-2189

# Spares/Consumables



#### Porous Plate 60mm Sq.

Code: 26-2181/10

#### **Alternatives**



### 2.42" Diameter Shear Box Assembly (California Type)

Code: 26-2223

+44 (0)1525 249 200 **= +1 (800) 323 1242** Page 2 of 3 www.ele.com

# **Product Sheet**







### **Shear Box Assembly 100mm Sq.**

Code: 26-2197



## Specimen Cutter 60mm Sq.

Code: 26-2185



## Specimen Extrusion Tool 60 mm sq.

Code: 26-2189

+44 (0)1525 249 200 **= +1 (800) 323 1242** www.ele.com

Page 3 of 3