Product Sheet







50MI Density Bottle with Perforated Stopper.

Code: 24-2900

Product Group: Particle Density: Pyknometer

Methods

Supplied complete with capillary vent stopper. Conforms to BS 733, and ISO 3507. Gay-Lussac

type.

Standards ASTM C128, AASHTO T84

Accessories



100MI Density Bottle with Perforated Stopper.

Code: 24-2950



Density Bottle 25Ml.

Code: 24-2890

Alternatives



100Ml Density Bottle with Perforated Stopper.

Code: 24-2950



Density Bottle 25Ml.

Code: 24-2890

+44 (0)1525 249 200

= +1 (800) 323 1242