



Heat of Hydration Apparatus

Code: [38-4600/01](#)

Product Group: [Heat of Hydration Cement](#)

Designed to determine the heat of hydration of cement to BS 4550.

Comprises of a dewar flask, internally lagged case, constant speed stirrer and filter funnel.

Please note thermometer must be purchased separately.

For 220 - 240 V AC, 50 Hz, 1 ph.

Standards

BS 4550

Further Information

****NOT AVAILABLE FOR SALE IN THE USA****

Accessories



200 mm Diameter Receiver

Code: [79-0020](#)



200mm Dia BS Sieve 125Mic Stainless Steel Mesh.

Code: [79-0110](#)



200mm Dia BS Sieve 600Mic Stainless Steel Mesh.

Code: [79-0200](#)



Desiccator - Non-vacuum

Code: [82-2100](#)



Digital Thermometer and Probe for 0.001 Dec C Resolution

Code: [38-4600/15](#)



Digital Thermometer and Probe for 0.01 Dec C Resolution

Code: [38-4600/16](#)



Electronic Balance - 210g x 0.0001g

Code: 78-5325/01



Electronic Balance - 220g x 0.001g

Code: 78-5335/01



Electronic Balance - 6 kg x 0.1g

Code: 78-5600



Harvard Trip Balance

Code: 78-7090



Sieve Brush Nylon

Code: 79-7210



Silica Gel

Code: [82-7091](#)



Spatula 100 mm

Code: [81-0100](#)



Weighing Bottle

Code: [82-1540](#)



Weights for Harvard Trip Balance

Code: [78-7110](#)

Spares/Consumables



Beckman Thermometer and Reader Glass

Code: [38-4600/12](#)



Dewar Flask

Code: [38-4600/10](#)



Glass Paddle for Stirrer

Code: [38-4600/13](#)

Alternatives



Mortar and Pestle Porcelain.

Code: [23-3505](#)



Muffle Furnace

Code: [83-4170/01](#)



Porcelain Crucible

Code: [82-3320](#)