

Sand Equivalent Value



Sand Equivalent Apparatus ASTM Aashto with Carrying Case.

Code: [24-4919](#)

Attributes: [Description](#) = Sand equivalent apparatus ASTM D2419

Standards: [ASTM D2419](#), [AASHTO T176](#)



5 Litre Syphon Assembly with Connection for 5mm Id Tubing.

Code: [24-4925](#)

Attributes: [Description](#) = Syphon assembly

Standards: [ASTM D2419](#), [AASHTO T176](#)



Calcium Chloride 2.5Kg.

Code: [24-4930](#)

Attributes: [Description](#) = Calcium chloride, [Weight](#), kg = 2.5

Standards: [ASTM D2419](#), [AASHTO T176](#)



Formaldehyde 40% Aqueous Solution 2.5 Litre.

Code: [24-4932](#)

Attributes: [Description](#) = Formaldehyde

Standards: [ASTM D2419](#), [AASHTO T176](#)



Glycerol Analar 2.5 Litre

Code: [24-4934](#)

Attributes: [Description](#) = Glycerol analar

Standards: [ASTM D2419](#), [AASHTO T176](#)



Motorised Sand Equivalent Shaker ASTM. 220-240VK 50Hz 1Ph.

Code: [24-4945/01](#)

Attributes: [Description](#) = Motorised sand equivalent shaker, [Power Supply](#) = 220 - 240 V AC, 50 Hz, 1 ph

Standards: