

Particle Density: Pyknometer Methods



50ml Density Bottle with Perforated Stopper.

Code: 24-2900

Attributes: Description = Density bottle, Capacity, ml = 50

Standards: ASTM C128, AASHTO T84



100ml Density Bottle with Perforated Stopper.

Code: 24-2950

Attributes: Description = Density bottle, Capacity, ml = 100

Standards: ASTM C128, AASHTO T84



Density Bottle 25ml.

Code: 24-2890

Attributes: Description = Density bottle, Capacity, ml = 25

Standards: EN 1097-7, ASTM D854, AASHTO T100, ISO 3507



Pyknometer Complete with Non-Corrodable Cone and Rubber Seal 1 kg Capacity

Code: 24-2885

Attributes: Description = Pyknometer, Gasket diameter, i.d; o.d " = 2-1/4" (57.2 mm) i.d.; 2-7/8" (73 mm) o.d.

Standards: ASTM C128, AASHTO T84
