

Seed technology




Stand for Jacobsen Seed Testing Bath.

Code: [542-102](#)

Attributes: [Description](#) = Stand for Jacobsen Seed Testing Bath

Standards:



Glass Plate for Jacobsen Seed Testing Bath. (13 Plates Required for Each Bath).

Code: [542-104](#)

Attributes: [Description](#) = Glass Plate for Jacobsen Seed Testing Bath

Standards:



Plastic Covers Pack of 10 for Seed Testing Bath (7 Packs Required for Each Bath).

Code: [542-106](#)

Attributes: [Description](#) = Plastic Covers Pack of 10 for Seed Testing Bath

Standards:



Germination Cabinet with sealed stainless steel chamber 220-240v 50hz 1ph

Code: [542-124/01](#)

Attributes: [Description](#) = Germination Cabinet

Standards:



Laboratory Seed Blower 220-240V 50Hz 1ph

Code: [542-044/01](#)

Attributes: [Description](#) = Laboratory seed blower

Standards:



Jacobsen Seed Testing Bath.

Code: [542-100/01](#)

Attributes: [Description](#) = ELE Jacobsen Seed Testing Bath

Standards:

Product Group Sheet

www.ele.com
+44 (0) 01525 249 200

