



IF IT'S WORTH BUILDING, IT'S WORTH TESTING

# Contents

[www.ele.com](http://www.ele.com)  
+44 (0) 01525 249 200



This wishlist was generated on 09/04/2019, and contains the following Products:

---

26-2275/01

**Laboratory Vane Apparatus (Motorised)**

---

26-2279

**Vane 12.7mm Diameter X 12.7mm Long.**

---

26-2281

**Vane 12.7mm Diameter X 19mm Long**

---

26-2283

**Vane 12.7mm Diameter X 25.4mm Long.**

---



## Laboratory Vane Apparatus (Motorised)

Code: [26-2275/01](#)

Product Group: [Soil Lathes](#), [Soil Trimmers \(Lathes\)](#)

Laboratory Vane Apparatus with 4 Springs and 12.7mm Vane; Motorised

Comes with 38 and 110mm sample holder

The equipment is based on a design by the Transport Research Laboratory, England, and is available in hand or motorised versions. Stress is applied through the 12.7 x 12.7 mm vane by means of any one of four calibrated springs. The motorised version produces a shearing rate of 10° per minute.

Standards Compliance: BS 1377-7 EN DD  
ENV1997-2



## Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)

Product Group: [Laboratory Vane Apparatus](#),  
[Laboratory Vane Apparatus](#)

## Standards

BS 1377

## Accessories



### Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



### Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

## Spares/Consumables



### Set of 4 Calibrated Springs for Laboratory Vane Apparatus.

Code: [26-2275/10](#)

## Alternatives



### Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



### Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)



## Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)

Product Group: [Laboratory Vane Apparatus](#),  
[Laboratory Vane Apparatus](#)

## Standards

BS 1377

## Accessories




### Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



### Vane 12.7mm Diameter X 19mm Long


Code: [26-2281](#)



### Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

## Spares/Consumables




## Set of 4 Calibrated Springs for Laboratory Vane Apparatus.

Code: [26-2275/10](#)

---

### Alternatives




## Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



## Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



## Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)



## Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

Product Group: [Laboratory Vane Apparatus](#),  
[Laboratory Vane Apparatus](#)

## Standards

BS 1377

## Accessories



### Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



### Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)




### Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

## Spares/Consumables






## Set of 4 Calibrated Springs for Laboratory Vane Apparatus.

Code: [26-2275/10](#)


---

### Alternatives



## Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



## Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



## Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)