



IF IT'S WORTH BUILDING, IT'S WORTH TESTING

Contents

www.ele.com
+44 (0) 01525 249 200



This wishlist was generated on 26/04/2019, and contains the following Products:

26-2270

Laboratory Vane Apparatus with 4 Springs and 12.7mm Vane; Hand Operated.



Laboratory Vane Apparatus with 4 Springs and 12.7mm Vane; Hand Operated.

Code: [26-2270](#)

Product Group: [Laboratory Vane Apparatus](#),
[Laboratory Vane Apparatus](#)

The equipment is based on a design by the Transport Research Laboratory, England, and is available in hand or motorised versions. Stress is applied through the 12.7 x 12.7 mm vane by means of any one of four calibrated springs.

Standards

BS 1377

Further Information

Hand-operated, complete with 4 calibrated springs and a 12.7 mm diameter x 12.7 mm long vane.

Accessories




Vane 12.7mm Diameter X 12.7mm Long.

Code: [26-2279](#)



Vane 12.7mm Diameter X 19mm Long

Code: [26-2281](#)



Vane 12.7mm Diameter X 25.4mm Long.

Code: [26-2283](#)

Spares/Consumables



Set of 4 Calibrated Springs for Laboratory Vane Apparatus.

Code: [26-2275/10](#)

Alternatives



Digital Direct/Residual Shear Apparatus Complete with Lever Loading Assembly. 110-120V 50/60Hz 1Ph.

Code: [26-2114/02](#)



Digital Direct/Residual Shear Apparatus Complete with Lever Loading Assembly. 220-240V 50/60Hz 1Ph.

Code: [26-2114/01](#)