



Pyknometer Complete with Non-Corrodable Cone and Rubber Seal 1 kg Capacity

Code: [24-2885](#)

Product Group: [Particle Density: Pyknometer
Methods, Pycnometer Top and Jar](#)

Glass jar complete with non-corrodable cone.

Standards

ASTM C128, AASHTO T84

Accessories



100MI Density Bottle with Perforated Stopper.

Code: [24-2950](#)



50MI Density Bottle with Perforated Stopper.

Code: [24-2900](#)



Density Bottle 25MI.

Code: [24-2890](#)



Pyknometer Gasket

Code: [24-2885/11](#)



Pyknometer Jar Only

Code: [24-2885/10](#)



Pyknometer Top

Code: [24-2885/12](#)

Alternatives



Sand Absorption Cone

Code: [42-1700](#)